

**NEW** | **HEINE K180 LED F.O. otoscope**  
**Standard instrument for a smaller budget**



Otoscope with fiber optic technology.  
Housing made of lightweight, shockproof polycarbonate.

**DETAILED VIEW**

- ▶ High-resolution and distortion-free images due to acrylic viewing window
- ▶ 3x magnification
- ▶ True-to-life colour rendering
- ▶ Even and reflex-free illumination due to fiber optic technology
- ▶ Particularly bright due to LED<sup>HQ</sup>

**ENHANCED DURABILITY**

- ▶ Sturdy polycarbonate housing
- ▶ Lightweight yet durable

**INTUITIVE OPERABILITY**

- ▶ Practical one-finger operation
- ▶ Always the ideal light intensity thanks to continuous brightness control
- ▶ Easy instrumentation due to fixed, swivelling viewing window
- ▶ Safe, hygienic reprocessing
- ▶ Pneumatic eardrum testing by means of an integrated insufflation port

**FLEXIBLE POWER SOURCE OPTIONS**

- ▶ Rechargeable battery: with USB or table charger
- ▶ EN 200 | EN 200-1 wall transformer
- ▶ Replaceable batteries
- ▶ **A tip for a better future.** HEINE AllSpec disposable tips – made of recycled plastic. Read more on [page 15](#)