

MEASURING SYSTEMS
AND SCALES SINCE 1840

seca
Precision for health

PRECISION
FOR
HEALTH

2024

seca full scale

Choosing seca means more than just selecting a scale.

It means choosing excellence, tradition, precision, innovation, reliability, security, inventiveness, comfort, and efficiency.

Choosing seca is choosing the future of medical measuring and weighing.



As one of the first Med Tech companies in Europe, seca has passed the official MDR (Medical Device Regulation (EU) 2017/745) audit even before the new EU regulation came into effect. All products developed by seca will be MDR certified within 2024.



1840

the birth year of the future world market leader: two decades ahead of the telephone's debut.

100%

progress powered by precision.

17

branches worldwide: we are your next-door tech wizards.

70+

in-house engineers and developers keep us ahead in innovation.

68%

market share: serving health care professionals around the globe.

175k

newborns touch a seca scale for their first time every day.

110+

exclusive partnerships: the VIP pass to the seca experience.



1st

scale that knows your name: we redefine personalization with ID workflows and automatic measurements.

EMR

validation: direct and secure data transmission. No errors. No fuss.

All seca products shown together with this logo are capable of transmitting measurements into an Electronic Medical Record (EMR) system.



180+

years of medical device expertise: our pioneer spirit still thrives.



24 Million

daily measurements powered on a seca device worldwide: seca's precision in action.

97%

accurate in muscle mass measurement compared to a whole-body MRI. Beat that!



SERVICE

our award-winning commitment to medical care.



2

CEOs, brothers and owners steering the family legacy into its fifth generation.



Product Overview

1 — Baby scales

2 — Pediatric measuring systems

3 — Measuring stations

4 — Column scales

5 — Flat scales

6 — Height and circumference measuring systems

7 — Handrail scales


8 — Wheelchair scales

9 — Bed and stretcher scales


10 — Chair scales

11 — Specialized scales

12 — Accessories



This is an interactive document. Click on the relevant product category to access it directly. You can get back to the Overview by clicking the "Back to product overview" button.



Baby scales

Babies who feel safe and secure are easier to weigh.

Baby scales by seca make weighing fast and comfortable, thanks to the skin-friendly materials that immediately take on and radiate body warmth. With designs that make the baby feel safe and advanced functions that shorten the weighing process, such as seca integration solutions. They allow you to transmit measurements quickly, securely and wirelessly to any EMR (Electronic Medical Record) system.

seca 727

EMR-validated baby scale with very precise graduation

- + Utmost accuracy: suitable for premature infants
- + Weighing tray made of heat-retaining plastic
- + Backlit, easy-to-read LCD display
- + Included rechargeable batteries allow mobile use without searching for a power outlet
- + Three-stage adjustable damping for fast weighing
- + EMR-validated: transmit measured data directly to an electronic medical record system or printer

Capacity	44 lbs / 20 kg
Graduation	0.1 oz / 1 g < 6 kg > 2 g
Power supply	Power adapter, Rechargeable battery
Interfaces	Wi-Fi / Ethernet optional with seca 452
Accessories	<ul style="list-style-type: none"> + Mobile cart seca 402 + Mobile cart with drawer seca 403 + Interface module seca 452 + Rechargeable battery seca 454



[Go to product website](#)

[Download Product Sheet](#)

[Back to product overview](#)

seca 333 i

EMR-validated baby scale with Wi-Fi function

- + EMR-validated via Wi-Fi: transmit measured data directly to an electronic medical record (EMR) system or printer
- + USB port for barcode scanner for user and patient identification
- + Optional with digital measuring rod seca 234 or analog measuring rod seca 232 n
- + Damping system for fast weighing
- + Easy-to-read, backlit LCD display

Capacity	44 lbs / 20 kg
Graduation	0.2 oz / 5 g
Power supply	Batteries, Power adapter
Interfaces	USB for barcode scanner, Wi-Fi
Accessories	<ul style="list-style-type: none"> + Analog measuring rod seca 232 n + Digital measuring rod seca 234 + Mobile cart seca 402 + Mobile cart with drawer seca 403 + Carrying case seca 428 + Rechargeable battery seca 454



seca 374

Baby scale with extra large weighing tray


- + Weighing tray with high side walls for increased safety
- + Weighing and measuring in one step thanks to optional length measuring rod
- + Three-step damping for fast weighing
- + Large tray is suitable for weighing toddlers while sitting
- + Medical-grade seca scale technology with precise graduation



Image contains optional measuring rod seca 233.



Capacity	44 lbs / 20 kg
Graduation	0.2 oz < 22 lbs > 0.5 oz / 5 g < 10 kg > 10 g
Power supply	Batteries, Power adapter (optional)
Accessories	<ul style="list-style-type: none"> + Analog measuring rod seca 233 + Switch-mode power adapter seca 400 + Mobile cart seca 402 + Mobile cart with drawer seca 403



[Go to product website](#)

[Download Product Sheet](#)

[Back to product overview](#)

seca 334


Mobile baby scale with flat weighing tray

- + Lightweight with a carrying handle for mobile use
- + Measuring and weighing can be done in a single step thanks to an optional length measuring rod
- + Flat construction for easy storage in drawers
- + Large, easy to read LCD display
- + Wide range of functions: TARE, BMIF and Auto-HOLD



Image contains optional measuring rod seca 232.

Capacity	44 lbs / 20 kg
Graduation	0.2 oz < 22 lbs > 0.5 oz / 5 g < 10 kg > 10 g
Power supply	Batteries, Power adapter (optional)
Accessories	<ul style="list-style-type: none"> + Analog measuring rod seca 232 + Switch-mode power adapter seca 400 + Mobile cart seca 402 + Mobile cart with drawer seca 403 + Carrying case seca 428



seca 354

2-in-1 mobile baby scale and flat scale for toddlers

- + Weigh babies lying down and toddlers standing up
- + Compact, robust and lightweight: ideal for mobile use
- + High-quality baby scale at an economical price
- + Wide range of functions: TARE, HOLD and Breast-Milk-Intake-Function (BMIF)



Capacity	44 lbs / 20 kg
Graduation	0.2 oz < 22 lbs > 0.5 oz / 5 g < 10 kg > 10 g
Power supply	Batteries
Accessories	<ul style="list-style-type: none"> + Mobile cart seca 402 + Mobile cart with drawer seca 403 + Backpack seca 409 + Carrying case seca 413 + Carrying case seca 414



[Go to product website](#)


[Download Product Sheet](#)

[Back to product overview](#)

seca 725

Mechanical baby scale with sliding weights

- + Classic construction with high-quality precision mechanics
- + Same level of accuracy as a digital baby scale
- + Corrosion resistant inside and outside: Moisture is not a problem
- + Removable weighing tray is easy to clean and perfect for mobile use
- + Additional sliding weight for the deduction of the tare weight of cloths, etc.
- + Depending on the model, with level and leveling feet



Capacity	32 lbs or 16 kg
Graduation	1/4 oz or 5 g

[Go to product website](#)

[Download Product Sheet](#)

[Back to product overview](#)

2

Pediatric measuring systems

Adaptable devices are preferred.

Clinics, medical practices and health care facilities have varied conditions and requirements. seca responds with a wide range of products for pediatrics – from infantometers to a portable measuring mat and from measuring rods to a head circumference measuring tape. seca has the right solution for every application, whether the devices are for in-patient or mobile use, in large or small spaces, for everyday or infrequent use. They all share the renowned seca quality which always lives up to the most demanding medical requirements.


seca 416

Infantometer for measuring babies and toddlers

- + Robust, extra large construction for precise measuring
- + Ensures correct lying position due to trough-shaped measuring surface
- + Measurement scale printed along top side for easy reading
- + Foot positioner cannot be moved by the infant thanks to the locking button



Measuring range	13 – 39 in / 33 – 100 cm
Graduation	1/16 in / 1 mm



[Go to product website](#)

[Download Product Sheet](#)

[Back to product overview](#)

seca 417

Light, space-saving and stable measuring board also ideal for mobile use

- + Simple operation for quick mobile use
- + Folds up for space-saving storage
- + High-quality materials for long service life
- + Smooth surfaces for easy cleaning



Measuring range	4 – 39 in / 10 – 100 cm
Graduation	1/16 in / 1 mm
Accessories	+ Backpack seca 409 + Carrying case seca 412



[Go to product website](#)


[Download Product Sheet](#)

[Back to product overview](#)

seca 210

Mobile measuring mat for babies and toddlers

- + Light-weight and easy to clean with most commercially available disinfectants
- + Extremely robust and durable
- + Easy-care, washable material, compatible with all common disinfectants
- + Measuring mat rolls up for easy transport and space-saving storage
- + Fixed head piece and the sliding foot positioner make it simple to use
- + Scale with cm and inches



Measuring range	4 – 39 in / 10 – 99 cm
Graduation	1/4 in / 5 mm
Accessories	+ Backpack seca 409



[Go to product website](#)

[Download Product Sheet](#)

[Back to product overview](#)


seca 207

Baby measuring rod with large calipers

- + Large head and foot positioner for fast and simple measuring
- + Secure table or wall mounting possible with screw fittings
- + Enormous space savings thanks to folding calipers
- + Made of aluminum for durability
- + Abrasion-resistant scale due to screen printing



Measuring range	3 – 39 in / 7 – 99 cm
Graduation	1/8 in / 1 mm



[Go to product website](#)


[Download Product Sheet](#)

[Back to product overview](#)

seca 212

Measuring tape for head circumference of babies and toddlers

- + Hygienic, disposable measuring tape of non-stretch Teflon material
- + Double function: head circumference with the front side of the tape and facial symmetry with the rear side
- + Easy-to-read, color-coded precise results
- + Practical dispenser (content: 15 tapes)



Measuring range 1 – 23 in / 3 – 59 cm

Graduation 1/16 in / 1 mm

[Go to product website](#)

[Download Product Sheet](#)

[Back to product overview](#)

3

Measuring stations

Measuring and weighing in a single step.

Measuring stations let you measure and weigh in one step to effectively determine patient height and weight. Be it an ultrasonic measuring station or measuring station with digital head slider, seca offers state-of-the-art technology, high-quality materials, elegant design and functions that make daily medical routine easier. EMR-validation means that the user no longer has to record the measured data manually and the seca Scale-up Line displays the name or ID of the patient.

seca Scale-up Line


EMR-validated ultrasonic measuring station with ID-Display™

- + Display shows user and patient names or IDs for increased patient safety
- + Outstanding precision and high-speed measuring thanks to perfected ultrasound technology
- + Large platform made of tested seca Bearclaw Glass® for absolute breakage resistance and perfect hygiene
- + Illuminated LED footprints in the weighing platform show patients the correct standing position for accurate height measurement
- + Voice guidance to support the measurement process
- + Convenient transport with extremely robust rubber wheels
- + Optional animal sticker facilitates measuring and weighing of children
- + EMR-validated: transmit measured data directly to an electronic medical record system or printer

↓ For technical data see next page



Image contains optional barcode scanner mount seca 463 and barcode scanner from a third-party vendor.



[Back to product overview](#)

seca Scale-up Line

Order code **ON SMM IUT NN**

Flat scale seca 654
+ Ultrasonic height measurement seca 257

Order code **ON SMM IVT NN**

Flat scale seca 654
+ Ultrasonic height measurement (short) seca 257 sh

Scale

Capacity	800 lbs / 360 kg
Graduation	0.1 lb / 50 g
Power supply	Power adapter
Interfaces	Ethernet, USB for barcode scanner, Wi-Fi
Accessories	+ Handrail seca 455 + Barcode scanner mount seca 463 + Animal sticker seca 487

Ultrasonic height measurement

Measuring range seca 257	1 ft 12 in – 7 ft 2 in / 60 – 220 cm
Measuring range seca 257 sh	1 ft 12 in – 6 ft 7 in / 60 – 200 cm
Graduation	1/8 in / 1 mm



[Go to product website](#)

[Download Product Sheet](#)

[Back to product overview](#)

seca 286

EMR-validated ultrasonic measuring station with voice guidance



- + Voice-guided self-measurement assists personnel
- + Outstanding precision and high-speed measuring thanks to perfected ultrasound technology
- + Hands-free measuring prevents the spread of germs
- + Precise 0.1 pounds graduation at a very high capacity of up to 660 pounds
- + Optional animal sticker facilitates measuring and weighing of children
- + EMR-validated: transmit measured data directly to an electronic medical record system or printer

Scale

Capacity	660 lbs / 300 kg
Graduation	0.1 lb / 50 g
Power supply	Power adapter
Interfaces	Wi-Fi / Ethernet optional with seca 452
Accessories	+ Interface module seca 452 + Animal sticker seca 487

Height measurement

Measuring range	1 ft 12 in – 6 ft 10 in / 60 – 210 cm
Graduation	1/8 in / 1 mm



[Go to product website](#)

[Download Product Sheet](#)

[Back to product overview](#)

seca 284

EMR-validated measuring station for height and weight

- + Measuring and weighing in one step with absolute precision
- + Precise 0.1 pound graduation at a very high capacity of up to 660 pounds
- + Precise head alignment by the seca Frankfurt Line Positioner
- + Three-line touchscreen for simultaneous weight, height and BMI display
- + High-quality glass platform with a heel positioner for exact three-point measurement
- + EMR-validated: transmit measured data directly to an electronic medical record system or printer

Scale

Capacity	660 lbs / 300 kg
Graduation	0.1 lb / 50 g
Power supply	Batteries for headpiece, Power adapter
Interfaces	Wi-Fi / Ethernet optional with seca 452
Accessories	+ Interface module seca 452 + Animal sticker seca 487

Ultrasonic height measurement

Measuring range	11 in – 7 ft 2 in / 30 – 220 cm
Graduation	1/8 in / 1 mm



[Go to product website](#)

[Download Product Sheet](#)

[Back to product overview](#)

4

Column scales

Indestructible support in patient care.

Thanks to their high precision, steadfast mechanics, durable aluminum, steel and cast iron construction as well as many sophisticated details, seca column scales are an indispensable part of the medical world. When used with a height measuring rod, they double as a measuring station and can determine weight, height and BMI in a single action. The column scales can be easily transported thanks to their transport wheels.

seca 703

EMR-validated column scale with 660 pounds capacity and measuring rod

- + Extremely robust, highly precise seca weighing technology with 660 pounds capacity
- + Non-tilting base made of heavy duty cast iron
- + Impact- and scratch-resistant powder coating
- + Can be used as a measuring station with included height measuring rod
- + Optional panda bear facilitates measuring and weighing of children
- + EMR-validated: transmit measured data directly to an electronic medical record system or printer

Scale

Capacity	660 lbs / 300 kg
Graduation	0.1 lb < 330 lbs > 0.2 lb / 50 g < 150 kg > 100 g
Power supply	Batteries, Power adapter (optional)
Interfaces	Wi-Fi / Ethernet optional with seca 452
Accessories	+ Switch-mode power adapter seca 400 + Interface module seca 452

Telescopic measuring rod seca 220 (included)

Measuring range	24 – 78 in / 60 – 200 cm
Graduation	1/8 in / 1 mm




Image contains included measuring rod seca 220.



[Go to product website](#)

[Download Product Sheet](#)

[Back to product overview](#)

seca 703 s

EMR-validated column scale with stadiometer and 660 pounds capacity


- + The height measuring rod with a heel plate is ideal for three-point measurement for patients specifically affected by obesity
- + Display faces the user for increased visual ergonomics
- + Extremely robust, highly precise seca weighing technology with 660 pounds capacity
- + Non-tilting base made of heavy duty cast iron
- + Impact- and scratch-resistant powder coating
- + EMR-validated: transmit measured data directly to an electronic medical record system or printer

Scale

Capacity	660 lbs / 300 kg
Graduation	0.1 lb < 330 lbs > 0.2 lb / 50 g < 150 kg > 100 g
Power supply	Batteries, Power adapter (optional)
Interfaces	Wi-Fi / Ethernet optional with seca 452
Accessories	+ Panda seca 459

Stadiometer (included)

Measuring range	3 – 90 in / 6 – 230 cm
Graduation	1/8 in / 1 mm



[Go to product website](#)

[Download Product Sheet](#)

[Back to product overview](#)

seca 797

EMR-validated column scale with eye-level display and Wi-Fi function


- + Medical-grade, robust scale for years of use with high patient volume
- + Digital measuring rod with an extra long measuring range
- + Two-line display shows height and weight
- + Damping system and Auto-HOLD guarantee quick and precise weighing
- + 8-year warranty (4x the industry standard)
- + Optional: Panda bear facilitates measuring and weighing of children
- + EMR-validated with built-in Wi-Fi module

Scale

Capacity	550 lbs / 250 kg
Graduation	0.2 lb / 100 g
Power supply	Batteries (optional), Power adapter
Interfaces	USB for barcode scanner, Wi-Fi
Accessories	+ Panda seca 459 + Handrail seca 477

Telescopic measuring rod (included)

Measuring range (total)	4 – 90 in / 10 – 230 cm
Measuring range (digital)	48 – 90 in / 122 – 230 cm
Graduation	1/8 in / 1 mm



[Go to product website](#)

[Download Product Sheet](#)

[Back to product overview](#)

seca 787

EMR-validated column scale with eye-level display


- + Medical-grade, robust scale for years of use with high patient volume
- + Digital measuring rod with an extra long measuring range
- + Two-line display shows height and weight
- + Damping system and Auto-HOLD guarantee quick and precise weighing
- + 8-year warranty (4x the industry standard)
- + Optional: Panda bear facilitates measuring and weighing of children
- + EMR validated: Send results directly into the EMR

Scale

Capacity	550 lbs / 250 kg
Graduation	0.2 lb / 100 g
Power supply	Batteries, Power adapter (optional)
Interfaces	Wi-Fi / Ethernet optional with seca 452
Accessories	<ul style="list-style-type: none"> + Power adapter 401 + Interface module seca 452 + Panda seca 459 + Handrail seca 477

Telescopic measuring rod (included)

Measuring range (total)	4 – 90 in / 10 – 230 cm
Measuring range (digital)	48 – 90 in / 122 – 230 cm
Graduation	1/8 in / 1 mm



seca 777

Digital column scale with eye-level display


- + Medical-grade, robust scale for years of use with high patient volume
- + Built-in measuring rod with 4 – 90 inch measuring range
- + USB interface for PC connectivity
- + Large rubber coated castors for excellent movability
- + Oversized, low-profile full metal platform
- + 8-year warranty (4x the industry standard)
- + Optional: handrail seca 477 provides extra patient security
- + Optional: panda bear facilitates measuring and weighing of children

Scale

Capacity	550 lbs / 250 kg
Graduation	0.2 lb / 100 g
Power supply	Batteries, Power adapter (optional)
Interfaces	USB
Accessories	+ Power adapter 401 + Panda seca 459 + Handrail seca 477

Telescopic measuring rod (included)

Measuring range	4 – 90 in / 10 – 230 cm
Graduation	1/8 in / 1 mm



[Go to product website](#)

[Download Product Sheet](#)

[Back to product overview](#)

seca 769

Digital column scale with BMI function



- + Particularly robust and compact column scale
- + Includes BMI function, TARE and HOLD
- + Robust weighing platform made of steel
- + Easily movable with transport wheels
- + Can be powered with batteries for mobile use or a power adapter for stationary use
- + Expandable to measuring station with an optional height measuring rod

Scale

Capacity	450 lbs / 200 kg
Graduation	0.2 lb / 100 g
Power supply	Batteries, Power adapter (optional)
Accessories	+ Telescopic measuring rod seca 220 + Switch-mode power adapter seca 400

Telescopic measuring rod seca 220 (optional)

Measuring range	24 – 78 in / 60 – 200 cm
Graduation	1/8 in / 1 mm



[Go to product website](#)

[Download Product Sheet](#)


[Back to product overview](#)

seca 719

Column scale with four weight sensors for precise results

- + Display on a high column for easy reading
- + Comfortable switching on by stepping on the platform
- + High stability for demanding users
- + More than 180 years of experience for maximum precision

Capacity	400 lbs / 180 kg
Graduation	0.2 lb / 100 g
Power supply	Batteries



[Go to product website](#)

[Download Product Sheet](#)

[Back to product overview](#)

seca 700

Mechanical column scale with eye-level beam

- + Easy reading and operation at eye level
- + Requires no batteries or other power supply
- + Same level of accuracy as a digital column scale
- + Non-tilting and highly robust base made of cast iron
- + Easy movable due to transport castors
- + Expandable to measuring station with an optional height measuring rod

Scale

Capacity	400 lbs / 160 kg
Graduation	0.25 lb / 100 g

Telescopic measuring rod seca 220 (optional)

Measuring range	24 – 78 in / 60 – 200 cm
Graduation	1/8 in / 1 mm




Image contains optional measuring rod seca 220.



[Go to product website](#)

[Download Product Sheet](#)

[Back to product overview](#)

seca 755

Mechanical column scale with BMI display and evaluation



- + Fast and easy BMI determination
- + Simultaneous evaluation of nutritional status
- + Clear circular scale and low platform
- + Robust steel construction, high-quality and durable
- + Corrosion resistant inside and outside: moisture is not a problem
- + Non-slip tread surface thanks to rubber pad
- + Always ready for use without power connection or batteries

Scale

Capacity	350 lbs / 160 kg
Graduation	1 lb / 500 g

Telescopic measuring rod seca 224 (included)

Measuring range	60 – 200 cm
Graduation	1 mm



[Go to product website](#)

[Download Product Sheet](#)

[Back to product overview](#)

5

Flat scales

Saving on space does not mean saving on patient benefits.

Where space is limited, there remains plenty of space for seca flat scales. Their wide array of functions make the day-to-day routines in medical care easier. Additional features include high weight capacity of up to 360 kg and long-battery life. Useful accessories such as carrying cases and a stand for a separate display make using seca flat scales even more enjoyable.

seca 650

EMR-validated flat scale with ID-Display™ and stable glass platform

- + Display shows names or IDs of users and patients for enhanced safety
- + Maximum flexibility in display positioning with universal mount for tabletop or wall mounting
- + Large platform made of tested seca Bearclaw Glass® for absolute breakage resistance and perfect hygiene
- + Glass surface is easy to clean and disinfect (more hygienic than plastic)
- + Adjustable damping for fast and precise weighing
- + EMR-validated: transmit measured data directly to an electronic medical record system or printer



Capacity	800 lbs / 360 kg
Graduation	0.1 lb / 50 g
Power supply	Power adapter
Interfaces	Ethernet, USB for barcode scanner, Wi-Fi



[Go to product website](#)

[Download Product Sheet](#)


[Back to product overview](#)

seca 634

EMR-validated flat scale with extra-large platform

- + Extra-large platform with 800 pounds capacity
- + Cable remote display for wall mounting or table installation
- + Includes BMI calculator and three-stage damping
- + Wide range of functions: TARE, pre-TARE, HOLD and Auto-HOLD
- + EMR-validated: transmit measured data directly to an electronic medical record system or printer

Capacity	800 lbs / 360 kg
Graduation	0.1 lb / 50 g
Power supply	Batteries, Power adapter (optional)
Interfaces	Wi-Fi / Ethernet optional with seca 452
Accessories	+ Switch-mode power adapter seca 400 + Interface module seca 452 + Stand for display element seca 472



[Go to product website](#)

[Download Product Sheet](#)

[Back to product overview](#)

seca 869

Flat scale with cable remote display

- + Cable remote display for wall mounting or table installation
- + Reliable and durable: seca quality that pays off
- + Leveling feet for strong stability
- + Can be powered with batteries for mobile use or a power adapter for stationary use
- + With TARE function, HOLD function and BMI computer

Capacity	550 lbs / 250 kg
Graduation	0.2 lb < 330 lbs > 0.5 lb / 100 g < 150 kg > 200 g
Power supply	Batteries, Power adapter (optional)
Accessories	<ul style="list-style-type: none"> + Switch-mode power adapter seca 400 + Backpack seca 409 + Carrying case seca 414 + Carrying case seca 421 + Adapter element seca 437



[Go to product website](#)

[Download Product Sheet](#)

[Back to product overview](#)

seca 874 +

Flat scale with foot switch and optional Bluetooth® interface

- + Developed for mobile use and protected against environmental influences with an IP53 protection code
- + Platform with an anti-slip surface and antimicrobial plastic for optimal hygiene
- + Extremely stable steel base with a heavy-duty weight capacity
- + Convenient and easy to operate with a single foot switch for efficient use
- + LC-display can be oriented towards the patient or doctor, offering convenience and flexibility
- + Mother+child function simplifies the weighing of children in the arms of the mother or father
- + Optional Bluetooth® interface: enables error-free data transmission

Capacity	440 lbs / 200 kg
Graduation	0.2 lb / 100 g
Power supply	Batteries, Power adapter (optional)
Interfaces	none
Accessories	+ Switch-mode power adapter seca 400 + Power adapter 401




seca 874

Flat scale with foot switches and double display

- + Double display lets the doctor and patient read simultaneously
- + Reliable and durable: seca quality that pays off
- + Can be operated in a standing position thanks to the foot switches
- + Can weigh 10,000 times with a single set of batteries
- + Mother+child function to weigh toddlers

Capacity	440 lbs / 200 kg
Graduation	0.1 lb < 330 lbs > 0.2 lb / 50 g < 150 kg > 100 g
Power supply	Batteries
Accessories	<ul style="list-style-type: none"> + Backpack seca 409 + Carrying case seca 414 + Carrying case seca 415 + Adapter element seca 437



[Go to product website](#)

[Download Product Sheet](#)


[Back to product overview](#)

seca 874 dr

Flat scale with customizable label

- + Customizable label builds trust
- + Cost-effective, reliable and durable: seca quality that pays off
- + Double display allows the doctor and patient to read simultaneously
- + Leveling feet for strong stability
- + Mother+child function for hassle-free weighing of toddlers

Capacity	440 lbs / 200 kg
Graduation	0.1 lb < 330 lbs > 0.2 lb / 50 g < 150 kg > 100 g
Power supply	Batteries
Accessories	<ul style="list-style-type: none"> + Switch-mode power adapter seca 400 + Backpack seca 409 + Carrying case seca 414 + Carrying case seca 415 + Adapter element seca 437



[Go to product website](#)

[Download Product Sheet](#)


[Back to product overview](#)

seca 876

Flat scale with foot switch

- + Precise weighing technology for medical use
- + Mother+child function to weigh toddlers while cradled on the mother's or father's arm
- + Can be operated in a standing position thanks to foot switch and tap-on function
- + Leveling feet for strong stability
- + Can weigh up to 10,000 times with a single set of batteries

Capacity	550 lbs / 250 kg
Graduation	0.2 lb < 330 lbs > 0.5 lb / 100 g < 150 kg > 200 g
Power supply	Batteries
Accessories	<ul style="list-style-type: none"> + Backpack seca 409 + Carrying case seca 414 + Carrying case seca 415 + Adapter element seca 437



[Go to product website](#)

[Download Product Sheet](#)

[Back to product overview](#)


seca 813 bt

Flat scale with Bluetooth® interface and high capacity

- + Solid steel frame, up to 440 lbs capacity
- + Non-slip rubber base for safe standing during measuring
- + Large LCD digits for easy reading
- + Transmits weight via Bluetooth®
- + Compatible with Google Fit and Apple Health for simple weight trend tracking



Capacity	440 lbs / 200 kg
Graduation	0.2 lb / 100 g
Power supply	Batteries
Interfaces	Bluetooth®



[Go to product website](#)

[Download Product Sheet](#)

[Back to product overview](#)


seca 813

Flat scale with high capacity

- + Solid steel frame, up to 440 lbs capacity
- + Non-slip rubber base for safe standing during measuring
- + Large LCD digits for easy reading
- + Step-off function for weighing without switching on the scale
- + More than 180 years of experience for maximum precision
- + From the global market leader for medical weighing and measuring



Capacity	440 lbs / 200 kg
Graduation	0.2 lb / 100 g
Power supply	Batteries



[Go to product website](#)

[Download Product Sheet](#)

[Back to product overview](#)


seca 803

Digital flat scale with high-quality two-component rubber surface

- + Two-component rubber coating offers higher security and comfort while standing on the scale
- + Four stable rubber feet provide a stable feeling
- + Large, easy-to-read 1.1 inch LCD digits
- + Optionally available in black or white
- + More than 180 years of experience for maximum precision
- + From the global market leader for medical weighing and measuring



Capacity	330 lbs / 150 kg
Graduation	0.2 lb / 100 g
Power supply	Batteries



[Go to product website](#)

[Download Product Sheet](#)

[Back to product overview](#)


seca 807

Digital personal scale with extra-flat dimensions

- + Four-weighing-cell technology for extra flat construction
- + Step-off function for weighing without needing to switch on the scale
- + Safety glass construction offers high impact and shock resistance
- + Large, easy-to-read 1.1 inch LCD digits
- + More than 180 years of experience for maximum precision
- + From the global market leader for medical weighing and measuring



Capacity	330 lbs / 150 kg
Graduation	0.2 lb / 100 g
Power supply	Batteries



[Go to product website](#)


[Download Product Sheet](#)

[Back to product overview](#)


seca 762

Mechanical scale with large dial

- + Robust, powder-coated steel housing provides an attractive and longlasting finish
- + Corrosion resistant inside and outside: Moisture is not a problem
- + Clip fixation system for easy rubber-mat exchange
- + Ready for use without electrical connection or battery
- + High-quality medical grade personal scale
- + A weighing experience that no digital scale can offer



Capacity	320 lbs / 150 kg
Graduation	1 lb / 500 g



[Go to product website](#)

[Download Product Sheet](#)


[Back to product overview](#)

seca 750

Mechanical scale with robust steel housing

- + Classic design with high-quality mechanical operation
- + Large, easy-to-read circular scale
- + Corrosion resistant inside and outside: moisture is not a problem
- + Anti-slip support for a secure position
- + Tested medical device quality
- + Rugged housing made of steel with powder-coated lacquer
- + Independent of batteries or power

Capacity	320 lbs / 150 kg
Graduation	1 lb / 1 kg


[Go to product website](#)
[Download Product Sheet](#)
[Back to product overview](#)

6

Height and circumference measuring systems

When the quality is right,
so are the measurements.



Measuring rods and circumference measuring tapes are relatively simple devices in routine medical work – however, they still have to fulfill challenging clinical requirements. That's why seca offers non-abrasive measurement systems and parts made of distortion-resistant anodized aluminum.

seca 274

Freestanding stadiometer with display on the headpiece

- + Free-standing stadiometer: can be set up anywhere in the room
- + Digital height measurement: comfortable and precise
- + Precise head alignment by the seca Frankfurt line positioner
- + High-quality glass platform with heel plate
- + Leveling feet provide a very stable base
- + Includes heel positioner for ideal three-point measurement

Measuring range	11 in – 7 ft 2 in / 30 – 220 cm
Graduation	1/8 in / 1 mm
Power supply	Batteries



[Go to product website](#)

[Download Product Sheet](#)



[Back to product overview](#)

seca 264

Wall-mounted stadiometer with display on the headpiece

- + Wall mounting: absolutely stable and wobble-free height measurement
- + Digital height measurement: comfortable and precise
- + Precise head alignment by the seca Frankfurt line positioner
- + Headpiece with a backlit digital display for easy reading
- + Non-slip rubber mat with heel positioner
- + Includes heel positioner for ideal three-point measurement

Measuring range	11 in – 7 ft 2 in / 30 – 220 cm
Graduation	1/8 in / 1 mm
Power supply	Batteries



[Go to product website](#)

[Download Product Sheet](#)



[Back to product overview](#)

seca 217

Portable stadiometer combinable with flat scale

- + Can be connected to a mobile measuring station with a flat scale
- + Wall spacer for additional stability
- + Robust platform for secure stance
- + Convenient reading thanks to the scale printed on the side and the reading window
- + The measuring rod can be dismantled and stored securely in the base for easy transport
- + Includes heel positioner for optimal three-point height measurement

Measuring range	8 in – 6 ft 9 in / 20 – 205 cm
Graduation	1/8 in / 1 mm
Accessories	<ul style="list-style-type: none"> + Backpack seca 409 + Carrying case seca 412 + Carrying case seca 414 + Adapter element seca 437



[Go to product website](#)

[Download Product Sheet](#)



[Back to product overview](#)

seca 213

Portable stadiometer

- + Simple and easy to set up – no wall fastening necessary
- + Large floor plate ensures stability
- + Includes heel positioner for ideal three-point measurement
- + Easy to read thanks to the scale printed on the side and the reading window
- + The measuring rod can be dismantled and stored securely in the base for easy transport
- + Wall spacer for additional stability

Measuring range	8 in – 6 ft 9 in / 20 – 205 cm
Graduation	1/8 in / 1 mm
Accessories	+ Carrying case seca 412



[Go to product website](#)

[Download Product Sheet](#)



[Back to product overview](#)

seca 222

Telescopic measuring rod with large measuring range

- + Designed specifically for wall mounting
- + With side print for greater ease of reading
- + Extra wide head gauge
- + Includes heel positioner for ideal three-point measurement
- + Durable aluminum measuring rod
- + Abrasion-resistant scale due to high-quality screen printing

Measuring range	2 – 90 in / 6 – 230 cm
Graduation	1/8 in / 1 mm



[Go to product website](#)

[Download Product Sheet](#)


[Back to product overview](#)

seca 216

Measuring rod for children and adults

- + Flexible measuring range options
- + Easy-to-read measurement scale
- + Adjustable head positioner
- + Choice of two units of measurement
- + Simple assembly

Measuring range	max. 55.3 in, selectable from 1.4 – 90 in / max. 138 cm, selectable from 3.5 – 230 cm
Graduation	1/16 in / 1 mm



[Go to product website](#)

[Download Product Sheet](#)

[Back to product overview](#)

seca 206

Rolling measuring tape for wall mounting

- + Space-saving, roll-up mechanism
- + Simple read-out of results through display window
- + Durable tape made of metal

Measuring range 0 – 87 in or 0 – 220 cm

Graduation 1/8 in or 1 mm



[Go to product website](#)

[Download Product Sheet](#)

[Back to product overview](#)

seca 201


Ergonomic circumference measuring tape

- + Retractable measuring tape for precise body measurements
- + Printed scale on side makes results easy-to-read
- + Convenient one-handed operation fasten start of tape to retraction mechanism
- + Compact and easy to store



Measuring range 0 – 80 in or 0 – 205 cm

Graduation 1/8 in or 1 mm



[Go to product website](#)

[Download Product Sheet](#)

[Back to product overview](#)

Handrail scales

Firm support for safe weighing.

Handrail scales from seca offer every patient a secure grip and are the best fall prevention when weighing. They are designed to be easily accessible and to weigh precisely, no matter where and how the patient is supported. Safe access and support options simplify workflows and save staff the hassle of careful positioning of the patient or adjustments to the scale before weighing. Practical features such as transport wheels and folding functions make the large scales transportable and stowable.

seca Scale-up Line

EMR-validated mobile handrail scale with digital measuring rod and ID-Display™

- + Designed for mobile use, record weight and height in patient rooms, at the bedside, or wherever required
- + Heavy-duty, robust design to withstand high-traffic environments
- + Integrated digital height rod for automatic readings on the display
- + Sturdy steel handrail to reduce fall risk and support nursing staff
- + Exceptional precision with 0.1 pound graduation, suitable for pediatric patients
- + Spacious and hygienic glass platform made of break-resistant seca Bearclaw Glass®
- + Customizable touchscreen display shows weight, height and patient/user IDs
- + EMR-validated, enabling seamless data transmission directly to an Electronic Medical Record (EMR) system or printer*

*EMR integration features are functional exclusively when the device is connected to the AC Adapter and are not applicable when powered by a mobile battery pack.

↓ For technical data see next page



NEW Optimized for mobility!

Handrail scales



seca Scale-up Line

Order code **ON SHM ITT NN**

Flat scale seca 654


+ Handrail seca 455

+ Digital telescopic measuring rod seca 254

Capacity	800 lbs / 360 kg
Graduation	0.1 lb / 50 g
Power supply	Power adapter
Interfaces	Ethernet, USB for barcode scanner, Wi-Fi
Accessories	+ Rechargeable battery pack seca 453 + Panda seca 459 + Barcode scanner mount seca 463

Digital telescopic measuring rod seca 254 (optional)

Measuring range (total)	4 – 90 in / 10 – 228 cm
Measuring range (digital)	48 – 90 in / 122 – 228 cm
Graduation	1/8 in / 1 mm



[Go to product website](#)

[Download Product Sheet](#)



[Back to product overview](#)

seca Scale-up Line

EMR-validated handrail scale with ID-Display™

- + Display shows user and patient names or IDs for increased patient safety
- + Data transmission controlled by touchscreen
- + Safe weighing with a stable handrail
- + Lightweight, sophisticated design and extremely robust thanks to steel handrail that is directly attached to the cast iron pedestal
- + Large platform made of tested seca Bearclaw Glass® for absolute breakage resistance and perfect hygiene
- + Convenient transport with extremely robust rubber wheels
- + Optional panda bear facilitates measuring and weighing of children
- + EMR-validated: transmit measured data directly to an electronic medical record system or printer

↓ For technical data see next page



[Back to product overview](#)


seca Scale-up Line

Order code **ON SHM INT NN**

Flat scale seca 654

+ Handrail seca 455

Capacity	800 lbs / 360 kg
Graduation	0.1 lb / 50 g
Power supply	Power adapter
Interfaces	Ethernet, USB for barcode scanner, Wi-Fi
Accessories	<ul style="list-style-type: none"> + Ultrasonic height measurement seca 257 + Ultrasonic height measurement (short) seca 257 sh + Panda seca 459 + Barcode scanner mount seca 463



[Go to product website](#)

[Download Product Sheet](#)

[Back to product overview](#)

seca Scale-up Line

EMR-validated handrail scale with ID-Display™ and height measurement

- + Display shows user and patient names or IDs for increased patient safety
- + Outstanding precision and high-speed measuring thanks to perfected ultrasound technology
- + Large platform made of tested seca Bearclaw Glass® for absolute breakage resistance and perfect hygiene
- + Illuminated LED footprints in the weighing platform show patients the correct standing position for accurate height measurement
- + Voice guidance to support the measurement process
- + Safe weighing with a stable handrail
- + Optional animal sticker facilitates measuring and weighing of children
- + EMR-validated: transmit measured data directly to an electronic medical record system or printer

↓ For technical data see next page



seca Scale-up Line

Order code **ON SHM IUT NN**

- Flat scale seca 654
- + Handrail seca 455
- + Ultrasonic height measurement seca 257

Order code **ON SHM IVT NN**

- Flat scale seca 654
- + Handrail seca 455
- + Ultrasonic height measurement (short) seca 257 sh

Capacity	800 lbs / 360 kg
Graduation	0.1 lb / 50 g
Power supply	Power adapter
Interfaces	Ethernet, USB for barcode scanner, Wi-Fi
Accessories	+ Barcode scanner mount seca 463 + Animal sticker seca 487

Ultrasonic height measurement

Measuring range seca 257	40 – 87 in / 100 – 220 cm
Measuring range seca 257 sh	40 – 79 in / 100 – 200 cm
Graduation	1/8 in / 1 mm



[Go to product website](#)

[Download Product Sheet](#)

[Back to product overview](#)

seca 644

EMR-validated handrail scale with generous platform

- + Safe weighing due to sturdy handrail
- + Large non-slip platform with 800 pounds capacity
- + 360° swivel display for easy reading
- + With a BMI calculator and damping system
- + EMR-validated: transmit measured data directly to an electronic medical record system or printer

Capacity	800 lbs / 360 kg
Graduation	0.1 lb / 50 g
Power supply	Power adapter, Rechargeable battery
Interfaces	Wi-Fi / Ethernet optional with seca 452
Accessories	+ Telescopic measuring rod seca 223 + Interface module seca 452

Telescopic measuring rod seca 223 (optional)

Measuring range	2 – 90 in / 6 – 230 cm
Graduation	1/8 in / 1 mm




Image contains optional measuring rod seca 223.



[Go to product website](#)

[Download Product Sheet](#)

[Back to product overview](#)

8

Wheelchair scales

Remove barriers with wheelchair scales from seca.

Weighing immobile patients can be a complicated process, but it doesn't have to be. The weighing solutions from seca offer the possibility of weighing patients in wheelchairs or even in hospital beds. The low-level, generously dimensioned platforms with ramps and other innovative functions guarantee comfortable, effortless and precise weighing.

seca XLine Base

EMR-validated wheelchair scale with optional third ramp

- + XXL dimensions and capacity to accommodate bariatric wheelchairs
- + Built to withstand the highest loads with a capacity of 1,100 lbs
- + Extra low profile ensures a smooth and safe transition for wheelchair access onto the platform
- + Third ramp available: upgrade your seca XLine and further improve accessibility
- + Tailor your experience by choosing between the basic display or the ID-Display™ with EMR-interface
- + Tare-by-Scanning: simplify the weighing process by taring out the weight of the wheelchair by scanning its barcode
- + Optional: handrail equipped with an antimicrobial coating to enhance safety or serve as a standing aid
- + EMR-validated, enabling seamless data transmission directly to an Electronic Medical Record (EMR) system or to a printer

↓ For technical data see next page



NEW

Wheelchair scales



seca XLine Base

Order codes	ON SWX INB NN	ON SWX INT NN
Platform scale (43 x 43 in) seca 604	•	•
+ Two ramps seca 404	•	•
+ Guard rail seca 446	•	•
+ Guard rail with castors seca 448	•	•
+ Display holder seca 476	•	•
+ Basic display seca 495	•	•
+ ID-Display™ seca 496		•



Capacity	1,100 lbs / 500 kg
Graduation	0.2 lb / 100 g
Power supply	Power adapter
Interfaces	Ethernet, USB for barcode scanner, Wi-Fi
Accessories	+ <u>Third ramp seca 405</u>



[Go to product website](#)

[Download Product Sheet](#)



[Back to product overview](#)

seca XLine Hold

EMR-validated wheelchair scale with handrail and optional third ramp

- + Handrail equipped with an antimicrobial coating to enhance safety or serve as a standing aid
- + XXL dimensions and capacity to accommodate bariatric wheelchairs
- + Built to withstand the highest loads with a capacity of 1,100 lbs
- + Extra low profile ensures a smooth and safe transition for wheelchair access onto the platform
- + Third ramp available: upgrade your seca XLine and further improve accessibility
- + Tailor your experience by choosing between the basic display or the ID-Display™ with EMR-interface
- + Tare-by-Scanning: simplify the weighing process by taring out the weight of the wheelchair by scanning its barcode
- + EMR-validated, enabling seamless data transmission directly to an Electronic Medical Record (EMR) system or to a printer

↓ For technical data see next page



NEW

Wheelchair scales



seca XLine Hold

Order codes	ON SWX IRB NN	ON SWX IRT NN
Platform scale (43 x 43 in) seca 604	•	•
+ Handrail seca 442	•	•
+ Two ramps seca 404	•	•
+ Guard rail seca 446	•	•
+ Basic display seca 495	•	•
+ ID-Display™ seca 496		•




Image contains optional third ramp seca 405.

Capacity	1,100 lbs / 500 kg
Graduation	0.2 lb / 100 g
Power supply	Power adapter
Interfaces	Ethernet, USB for barcode scanner, Wi-Fi
Accessories	+ Third ramp seca 405



[Go to product website](#)

[Download Product Sheet](#)

[Back to product overview](#)

seca 684

EMR-validated multifunctional scale

- + Safe weighing of standing, sitting or wheelchair patients
- + Firmly anchored steel handrail for secure weighing of frail patients
- + Large, easy to access platform with non-slip surface
- + Wheel barrier provides extra safety
- + 360° swivel display for easy reading
- + Mobile due to integrated transport wheels
- + EMR-validated: transmit measured data directly to an electronic medical record system or printer



Capacity	800 lbs / 360 kg
Graduation	0.1 lb / 50 g
Power supply	Power adapter
Interfaces	Wi-Fi / Ethernet optional with seca 452
Accessories	+ Interface module seca 452



[Go to product website](#)

[Download Product Sheet](#)

[Back to product overview](#)


seca 676

EMR-validated wheelchair scale with handrail

- + Large wheelchair accessible platform
- + Convenient and foldable to save space
- + Mobile due to transport wheels
- + Handrail as a standing aid
- + Non-slip rubber mat
- + EMR-validated: transmit measured data directly to an electronic medical record system or printer



Capacity	800 lbs / 360 kg
Graduation	0.1 lb / 50 g
Power supply	Power adapter
Interfaces	Wi-Fi / Ethernet optional with seca 452
Accessories	+ Interface module seca 452



[Go to product website](#)

[Download Product Sheet](#)

[Back to product overview](#)

seca 664

EMR-validated wheelchair scale

- + Extra-large, easily accessible platform: ideal for wheelchairs
- + High capacity of up to 800 pounds
- + Side rails on the platform ensure safe rolling of a wheelchair onto and from the scale
- + Easy to fold together and transport
- + Non-slip rubber mat
- + EMR-validated: transmit measured data directly to an electronic medical record system or printer

Capacity	800 lbs / 360 kg
Graduation	0.1 lb / 50 g
Power supply	Power adapter, Rechargeable battery
Interfaces	Wi-Fi / Ethernet optional with seca 452
Accessories	+ Interface module seca 452



[Go to product website](#)

[Download Product Sheet](#)



[Back to product overview](#)

seca 674

EMR-validated platform scale with remote display

- + Medical-grade, robust scale for years of use with high patient volume
- + Flexible cable remote display for wall mounting or use with a stand
- + Extra large platform: suitable for wheelchairs or patients with obesity
- + Permanently save the weight of wheelchairs or standing aids with the Pre-TARE function
- + Low profile platform with ramp for easy rolling onto the scale
- + EMR-validated: transmit measured data directly to an electronic medical record system or printer


Capacity	800 lbs / 360 kg
Graduation	0.1 lb / 50 g
Power supply	Power adapter
Interfaces	Wi-Fi / Ethernet optional with seca 452
Accessories	+ Interface module seca 452 + Stand for display element seca 472



[Go to product website](#)

[Download Product Sheet](#)

[Back to product overview](#)



Bed and stretcher scales

Precise results in all situations.

seca bed and stretcher scales enable you to obtain patient weight with utmost precision even in challenging scenarios such as emergency arrivals or wherever bed-bound patients have to be weighted. Simply roll the hospital bed, stretcher, or gurney onto the weighing platform, use the pre-TARE function, and determine the patient's weight.



seca XLine Stretch

EMR-validated scale for hospital beds and stretchers

- + XXL dimensions to accommodate hospital beds, stretchers and gurneys with ease
- + Built to withstand the highest loads with a capacity of 2,200 lbs
- + Weigh patients quickly and safely during emergency admissions or dialysis
- + Eliminate the need for estimating the weight of bedridden patients, weigh them comfortable and precise
- + Tare-by-Scanning: simplify the weighing process by taring out the weight of the stretcher or bed by scanning its barcode
- + Extra low profile ensures a smooth and safe transition for pushing beds and stretchers onto the platform
- + After identification with a barcode scanner the ID-Display™ shows patient's ID, user ID and DOB
- + EMR-validated, enabling seamless data transmission directly to an Electronic Medical Record (EMR) system or to a printer

Capacity	2,200 lbs / 1,000 kg
Graduation	0.2 lb / 100 g
Power supply	Power adapter
Interfaces	Ethernet, USB for barcode scanner, Wi-Fi


[Go to product website](#)
[Download Product Sheet](#)


Order code **ON SBX INT NN**

- 2x Platform scale (43 x 43 in) seca 604
- + Two ramps seca 404
- + 2x Guard rail seca 446
- + 2x Guard rail with castors seca 448
- + Connection kit seca 438
- + Display holder seca 476
- + ID-Display™ seca 496


[Back to product overview](#)


seca 656

EMR-validated platform scale for stretchers

- + Weigh patients quickly and safely while lying down on the stretcher
- + Stable construction with very high capacity
- + Transportable due to wheels and integrated handles
- + Pre-TARE function for determining precise weight of patient
- + EMR-validated: transmit measured data directly to an electronic medical record system or printer



Capacity	800 lbs / 360 kg
Graduation	0.1 lb / 50 g
Power supply	Power adapter
Interfaces	Wi-Fi / Ethernet optional with seca 452
Accessories	+ Interface module seca 452 + Stand for display element seca 472



[Go to product website](#)

[Download Product Sheet](#)

[Back to product overview](#)

10

Chair scales

Comfortable and save
weighing while seated.

How can you gently and carefully weigh bedridden patients or people who are unsteady on their feet? The best choice is a seca chair scale, which you can take directly to the patient. That saves you – and nursing staff – from unnecessary steps.

seca 954

EMR-validated chair scale with precise graduation

- + 0.1 pound graduation for exact measurement while sitting
- + Fold-up arm and footrests for easy access
- + Ergonomic seat for increased patient comfort
- + Suitable for patients affected by obesity with up to 660 pounds weighing capacity
- + Sophisticated engineering and robust construction for a long lifespan
- + Parking brake for secure position
- + EMR-validated: transmit measured data directly to an electronic medical record system or printer

Capacity	660 lbs / 300 kg
Graduation	0.1 lb < 330 lbs > 0.2 lb / 50 g < 150 kg > 100 g
Power supply	Power adapter, Rechargeable battery
Interfaces	Wi-Fi / Ethernet optional with seca 452
Accessories	+ Interface module seca 452 + Rechargeable battery seca 454



[Go to product website](#)

[Download Product Sheet](#)

[Back to product overview](#)


seca 952

Chair scale for weighing while seated

- + The simple and secure solution for patients with decreased mobility
- + Fold-up arm and footrests make sitting down easy
- + Ergonomic seat for increased patient comfort
- + Wheels include brakes to ensure the scale stays in place when needed
- + Highest quality for a long service life



Capacity	440 lbs / 200 kg
Graduation	0.2 lb / 100 g
Power supply	Batteries



[Go to product website](#)

[Download Product Sheet](#)

[Back to product overview](#)

A person wearing green scrubs and white gloves is using a scale. The scale is a flat, rectangular platform with a digital display and buttons on the left side. The person is holding a metal tray in their left hand and has their right hand on the scale's surface. The scale is placed on a white surface. In the background, there are two large red 'T' symbols on a white background with a grid of thin lines.

Specialized scales

Small size, big help.

As a specialist for medical measuring and weighing, seca naturally offers solutions for special applications. For instance, the splash-protected organ and diaper scale and the extremely precise portion and diet scale for clinical nutritionists.


seca 856

Organ and diaper scale with stainless steel cover

- + Complete stainless steel-surface reduces contamination
- + Easy to clean and disinfect
- + Scale completely sealed from moisture
- + Pre-TARE function for permanent use with the container of your choice
- + Increase in metal housing protects display from liquid



Capacity	5 kg
Graduation	1 g < 3 kg > 2 g
Power supply	Batteries



[Go to product website](#)

[Download Product Sheet](#)

[Back to product overview](#)

seca 852

Portion and diet scale

- + Gram-precise display for the accurate weighing of food
- + Scale can be stowed in compact bowl
- + Food-safe plastic bowl is easy to clean
- + High 6.6 lbs capacity



Capacity	6.6 lbs / 3 kg
Graduation	0.05 oz / 1 g
Power supply	Batteries



[Go to product website](#)

[Download Product Sheet](#)

[Back to product overview](#)

12

Accessories

Just the right enhancements.

Enhance your seca products with practical accessories. For example, choose from measuring rods, power adapters, matching transport bags and so much more.





seca 402

Mobile cart for seca baby scales

- + Custom-made indentations in the surface hold seca baby scales firmly in place
- + Ideal working height of 31.5 inches
- + Optionally available paper roll holder helps you keep the work area clean

Dimensions (WxHxD)	28.1 × 31.5 × 17.1 in / 715 × 800 × 435 mm
Net weight	25.4 lbs / 11.5 kg

[Go to product website](#)

seca 403

Mobile cart with drawer for seca baby scales

- + Custom-made indentations in the surface hold seca baby scales firmly in place
- + Ideal working height
- + Integrated drawer and shelf offer extra storage space so you always have everything within reach

Dimensions (WxHxD)	28.1 × 31.5 × 17.1 in / 715 × 800 × 435 mm
Net weight	60.6 lbs / 27.5 kg

[Go to product website](#)

seca 452

Interface module to transmit measurements to EMR systems

- + Data transmission via WLAN or LAN
- + Patient identification right on device
- + No transcription errors

Dimensions (WxHxD)	3.7 × 4.4 × 1.4 in / 95 × 113 × 35 mm
Net weight	0.4 lbs / 0.18 kg

[Go to product website](#)

[Back to product overview](#)



seca 234

Digital measuring rod for baby scales
seca 333i

- + Measuring results are easy to read because of the illuminated LCD display
- + Includes head and foot positioners
- + Durable measuring rod made of aluminum
- + Ergonomically shaped adjustment which prevents the rod from tilting

Dimensions (WxHxD)	34.5 × 3.6 × 11.7 in / 877 × 92 × 297 mm
Net weight	2.9 lbs / 1.3 kg
Measuring range	14 – 32 in / 35 – 80 cm
Graduation	1/8 in / 1 mm

[Go to product website](#)



seca 232 n

Analog measuring rod for baby scales
seca 333i

- + Measuring results are easy to read with enhanced color markings
- + Includes head and foot positioners
- + Durable measuring rod made of aluminum
- + Scratch-resistant height rod made with a high-quality screen printing process

Dimensions (WxHxD)	34.5 × 3.6 × 10.5 in / 877 × 92 × 266 mm
Net weight	2.6 lbs / 1.2 kg
Measuring range	14 – 32 in / 35 – 80 cm
Graduation	1/8 in / 1 mm

[Go to product website](#)



seca 233

Measuring rod for baby scale seca 374

- + Measuring results are easy to read with enhanced color markings
- + Includes head and foot positioners
- + Durable measuring rod made of aluminum
- + Scratch-resistant height rod made with a high-quality screen printing process

Dimensions (WxHxD)	25.2 × 7 × 11.7 in / 641 × 179 × 296 mm
Net weight	1.5 lbs / 0.7 kg
Measuring range	14 – 32 in / 35 – 80 cm
Graduation	1/8 in / 1 mm

[Go to product website](#)

[Back to product overview](#)



seca 232

Analog measuring rod for baby scale
seca 334

- + Measuring results are easy to read with enhanced color markings
- + Includes head and foot positioners
- + Durable measuring rod made of aluminum
- + Scratch-resistant height rod made with a high-quality screen printing process

Dimensions (WxHxD)	25.2 × 5.7 × 10.9 in / 640 × 145 × 277 mm
Net weight	1.1 lbs / 0.51 kg
Measuring range	14 – 32 in / 35 – 80 cm
Graduation	1/8 in / 1 mm

seca 400

Switch-mode power adapter for seca
scales

Dimensions (WxHxD)	1.7 × 1.4 × 2.9 in / 42 × 35 × 73 mm
Net weight	0.3 lbs / 0.16 kg
Cable length	72.8 in / 1.9 m

seca 454

Rechargeable battery for mobile use of
EMR-integrated seca scales

Net weight	6.7 lbs / 3.05 kg
Cable length	45.2 in / 1.2 m

[Go to product website](#)

[Go to product website](#)

[Back to product overview](#)



seca 428

Large carrying case for seca baby scales

- + Provides convenient and comfortable transport with adjustable shoulder strap
- + Convenient. Comfortable transport with adjustable shoulder strap
- + High-quality: made of robust nylon
- + Transport safety. Water and dust resistant

Dimensions (WxHxD)	25.6 × 12.6 × 4.3 in / 650 × 320 × 110 mm
Net weight	2.6 lbs / 1.2 kg

[Go to product website](#)



seca 409

Large backpack for portable seca measuring and weighing solutions

- + Versatile. Suitable for seca baby scales, flat scales or other measuring systems, e.g. measuring boards
- + Versatile. Use the seca 409 as bag or backpack for safe and easy transport
- + Comfortable. Transport made easy by two strong handles, safety belt and padded straps
- + High-quality. Made of durable material – waterproof and dustproof

Dimensions (WxHxD)	14.6 × 24.8 × 5.5 in / 370 × 630 × 140 mm
Net weight	3.5 lbs / 1.6 kg

[Go to product website](#)



seca 414

Large carrying case for portable seca measuring and weighing solutions

- + Convenient. Comfortable transport with adjustable shoulder strap
- + High-quality: made of robust nylon
- + Transport safety. Water and dust resistant

Dimensions (WxHxD)	24.6 × 13.8 × 8.3 in / 625 × 350 × 210 mm
Net weight	6.2 lbs / 2.8 kg

[Go to product website](#)

[Back to product overview](#)



seca 413

Carrying case for baby seca scales

- + Convenient. Comfortable transport with adjustable shoulder strap
- + High-quality: made of robust nylon
- + Transport safety. Water and dust resistant

Dimensions (WxHxD)	22.4 × 5.9 × 13 in / 570 × 150 × 330 mm
Net weight	1.9 lbs / 0.9 kg

[Go to product website](#)[Back to product overview](#)



seca 409

Large backpack for portable seca measuring and weighing solutions

- + Versatile. Suitable for seca baby scales, flat scales or other measuring systems, e.g. measuring boards
- + Versatile. Use the seca 409 as bag or backpack for safe and easy transport
- + Comfortable. Transport made easy by two strong handles, safety belt and padded straps
- + High-quality. Made of durable material – waterproof and dustproof

Dimensions (WxHxD)	14.6 × 24.8 × 5.5 in / 370 × 630 × 140 mm
Net weight	3.5 lbs / 1.6 kg

[Go to product website](#)



seca 414

Large carrying case for portable seca measuring and weighing solutions

- + Convenient. Comfortable transport with adjustable shoulder strap
- + High-quality: made of robust nylon
- + Transport safety. Water and dust resistant

Dimensions (WxHxD)	24.6 × 13.8 × 8.3 in / 625 × 350 × 210 mm
Net weight	6.2 lbs / 2.8 kg

[Go to product website](#)



seca 412

Carrying case for portable seca measuring and weighing solutions

- + Convenient. Comfortable transport with adjustable shoulder strap
- + High-quality: made of robust nylon
- + Transport safety. Water and dust resistant

Dimensions (WxHxD)	27 × 16.8 × 2.2 in / 686 × 427 × 56 mm
Net weight	1.4 lbs / 0.65 kg

[Go to product website](#)

[Back to product overview](#)



seca 455

Handrail for seca Scale-up Line

Dimensions (WxHxD)	31.5 × 44.1 × 21.3 in / 801 × 1,120 × 540 mm
Net weight	15.4 lbs / 7 kg

[Go to product website](#)



seca 463

Barcode scanner mount*

- + Ensures a barcode scanner is always ready at hand
- + Straightforward and quick mounting on the column
- + Suitable for all standard barcode scanners

*Barcode scanner from a third-party vendor

Dimensions (WxHxD)	3 × 3 × 6 in / 76 × 75 × 152 mm
Net weight	0.4 lbs / 0.2 kg

[Go to product website](#)



seca 487

Animal sticker for seca ultrasonic measuring stations

- + Colourful animals to make measuring stations more child-friendly
- + Animal figures at eye-level help obtain a correct head position for height measurement
- + High quality: sticker stays colorful for a long time
- + Approved medical device: Robust and resistant to disinfectants
- + Self-adhesive for easy application

Dimensions (WxHxD)	2.4 × 92.4 × 0 in / 60 × 2,346 × 0.5 mm
Net weight	0.4 lbs / 0.2 kg

[Go to product website](#)

[Back to product overview](#)



seca 452

Interface module to transmit measurements to EMR systems

- + Data transmission via WLAN or LAN
- + Patient identification right on device
- + No transcription errors

Dimensions (WxHxD)	3.7 × 4.4 × 1.4 in / 95 × 113 × 35 mm
--------------------	---------------------------------------

Net weight	0.4 lbs / 0.18 kg
------------	-------------------

[Go to product website](#)

[Back to product overview](#)



seca 220

Telescopic measuring rod for seca column scales

- + Measurements are easy to read, thanks to printing on side of rod and color-coded markings
- + Aluminum makes rod especially durable
- + Abrasion-resistant scale made with high-quality screen printing process

Dimensions (WxHxD)	2.2 × 58.3 × 11.3 in / 55 × 1,480 × 286 mm
Net weight	1.5 lbs / 0.7 kg
Measuring range	24 – 78 in / 60 – 200 cm
Graduation	1/8 in / 1 mm

[Go to product website](#)



seca 459

Panda bear Bao Bao for seca measuring stations and column scales

- + The cute little panda makes the measuring process more appealing and child friendly
- + Quick and easy installation
- + Easy to clean, with rounded edges
- + Long durability through high-quality DIBOND® material
- + Compatible with seca 703, seca mBCA 554, seca 777, seca 787, and seca 797

Dimensions (WxHxD)	14.9 × 26.3 × 1.7 in / 380 × 670 × 44 mm
Net weight	1.1 lbs / 0.5 kg

[Go to product website](#)



seca 477

Integrated handrail for seca 777 series scales

- + Enhanced patient safety without affecting measurement
- + Extra support and stability during the weighing process
- + Enables a firm stance for patients in need of assistance
- + Easy to mount, thanks to predrilled holes in the platform
- + Compatible on seca series 777, 787, 797

Dimensions (WxHxD)	16.9 × 44.2 × 12.2 in / 429 × 1,122 × 309 mm
Net weight	8.8 lbs / 4 kg

[Go to product website](#)

[Back to product overview](#)



seca 452

Interface module to transmit measurements to EMR systems

- + Data transmission via WLAN or LAN
- + Patient identification right on device
- + No transcription errors

Dimensions (WxHxD)	3.7 × 4.4 × 1.4 in / 95 × 113 × 35 mm
Net weight	0.4 lbs / 0.18 kg

[Go to product website](#)



seca 400

Switch-mode power adapter for seca scales

Dimensions (WxHxD)	1.7 × 1.4 × 2.9 in / 42 × 35 × 73 mm
Net weight	0.3 lbs / 0.16 kg
Cable length	72.8 in / 1.9 m

[Go to product website](#)



401

Switch-mode power adapter for EMR-integrated seca scales and measuring stations

Dimensions (WxHxD)	3.1 × 1.6 × 2.9 in / 80 × 40 × 75 mm
Net weight	0.2 lbs / 0.1 kg
Cable length	78.7 in / 2 m

[Back to product overview](#)



seca 472

Stand for display elements of seca scales

- + Stable platform with non-slip rubber feet
- + Sturdy, durable construction
- + Can also be used as a stand for external printer

Dimensions (WxHxD)	18.1 × 40.7 × 18.1 in / 460 × 1,035 × 460 mm
Net weight	17.2 lbs / 7.8 kg

[Go to product website](#)



seca 437

Adapter element for stadiometer seca 217 for connection with seca flat scales

- + Converts the seca 217 and a flat scale in a complete mobile measuring station
- + An integrated click system secures a stable and safety weighing process
- + Easy and fast assembly
- + High-quality: made of robust plastic

Dimensions (WxHxD)	13.1 × 3.4 × 12.2 in / 333 × 86 × 310 mm
Net weight	1.3 lbs / 0.6 kg

[Go to product website](#)



seca 217

Portable stadiometer combinable with flat scale

- + Can be connected to a mobile measuring station with a flat scale
- + Wall spacer for additional stability
- + Robust platform for secure stance
- + Includes heel positioner for optimal three-point height measurement

Dimensions (WxHxD)	12.9 × 84.4 × 22.6 in / 328 × 2,145 × 574 mm
Net weight	7.9 lbs / 3.6 kg
Measuring range	8 in – 6 ft 9 in / 20 – 205 cm
Graduation	1/8 in / 1 mm

[Go to product website](#)

[Back to product overview](#)



seca 409

Large backpack for portable seca measuring and weighing solutions

- + Versatile. Suitable for seca baby scales, flat scales or other measuring systems, e.g. measuring boards
- + Versatile. Use the seca 409 as bag or backpack for safe and easy transport
- + Comfortable. Transport made easy by two strong handles, safety belt and padded straps
- + High-quality. Made of durable material – waterproof and dustproof

Dimensions (WxHxD)	14.6 × 24.8 × 5.5 in / 370 × 630 × 140 mm
Net weight	3.5 lbs / 1.6 kg

[Go to product website](#)



seca 414

Large carrying case for portable seca measuring and weighing solutions

- + Convenient. Comfortable transport with adjustable shoulder strap
- + High-quality: made of robust nylon
- + Transport safety. Water and dust resistant

Dimensions (WxHxD)	24.6 × 13.8 × 8.3 in / 625 × 350 × 210 mm
Net weight	6.2 lbs / 2.8 kg

[Go to product website](#)



seca 421

Spacious carrying case for seca flat scales

- + Provides convenient and comfortable transport with adjustable shoulder strap
- + Convenient. Comfortable transport with adjustable shoulder strap
- + High-quality: made of robust nylon
- + Transport safety. Water and dust resistant

Dimensions (WxHxD)	16.9 × 13.8 × 5.1 in / 430 × 350 × 130 mm
Net weight	1.8 lbs / 0.81 kg

[Go to product website](#)

[Back to product overview](#)



seca 415

Carrying case for the flat scales

- + Convenient. Comfortable transport with adjustable shoulder strap
- + High-quality: made of robust nylon
- + Transport safety. Water and dust resistant

Dimensions (WxHxD)	13 × 15 × 2.8 in / 330 × 380 × 70 mm
Net weight	1.8 lbs / 0.8 kg

[Go to product website](#)



seca 400

Switch-mode power adapter for seca scales

Dimensions (WxHxD)	1.7 × 1.4 × 2.9 in / 42 × 35 × 73 mm
Net weight	0.3 lbs / 0.16 kg
Cable length	72.8 in / 1.9 m

[Go to product website](#)



401

Switch-mode power adapter for EMR-integrated seca scales and measuring stations

Dimensions (WxHxD)	3.1 × 1.6 × 2.9 in / 80 × 40 × 75 mm
Net weight	0.2 lbs / 0.1 kg
Cable length	78.7 in / 2 m

[Back to product overview](#)



seca 452

Interface module to transmit measurements to EMR systems

- + Data transmission via WLAN or LAN
- + Patient identification right on device
- + No transcription errors

Dimensions (WxHxD)	3.7 × 4.4 × 1.4 in / 95 × 113 × 35 mm
--------------------	---------------------------------------

Net weight	0.4 lbs / 0.18 kg
------------	-------------------

[Go to product website](#)[Back to product overview](#)



seca 409

Large backpack for portable seca measuring and weighing solutions

- + Versatile. Suitable for seca baby scales, flat scales or other measuring systems, e.g. measuring boards
- + Versatile. Use the seca 409 as bag or backpack for safe and easy transport
- + Comfortable. Transport made easy by two strong handles, safety belt and padded straps
- + High-quality. Made of durable material – waterproof and dustproof

Dimensions (WxHxD)	14.6 × 24.8 × 5.5 in / 370 × 630 × 140 mm
Net weight	3.5 lbs / 1.6 kg

[Go to product website](#)



seca 414

Large carrying case for portable seca measuring and weighing solutions

- + Convenient. Comfortable transport with adjustable shoulder strap
- + High-quality: made of robust nylon
- + Transport safety. Water and dust resistant

Dimensions (WxHxD)	24.6 × 13.8 × 8.3 in / 625 × 350 × 210 mm
Net weight	6.2 lbs / 2.8 kg

[Go to product website](#)



seca 412

Carrying case for portable seca measuring and weighing solutions

- + Convenient. Comfortable transport with adjustable shoulder strap
- + High-quality: made of robust nylon
- + Transport safety. Water and dust resistant

Dimensions (WxHxD)	27 × 16.8 × 2.2 in / 686 × 427 × 56 mm
Net weight	1.4 lbs / 0.65 kg

[Go to product website](#)

[Back to product overview](#)



seca 437

Adapter element for stadiometer
seca 217 for connection with seca flat
scales

- + Converts the seca 217 and a flat scale in a complete mobile measuring station
- + An integrated click system secures a stable and safety weighing process
- + Easy and fast assembly
- + High-quality: made of robust plastic

Dimensions (WxHxD)	13.1 × 3.4 × 12.2 in / 333 × 86 × 310 mm
Net weight	1.3 lbs / 0.6 kg

[Go to product website](#)

[Back to product overview](#)



seca 257

Ultrasonic height measurement for seca Scale-up Line and seca mBCA 555/554

Dimensions (WxHxD)	6.8 × 90.7 × 17.1 in / 173 × 2,304 × 435 mm
Net weight	14.3 lbs / 6.5 kg
Measuring range	1 ft 12 in – 7 ft 2 in / 60 – 220 cm
Graduation	1/8 in / 1 mm

[Go to product website](#)



seca 257 sh

Ultrasonic height measurement (short) for seca Scale-up Line and seca mBCA 555/554

Dimensions (WxHxD)	6.8 × 82.8 × 17.1 in / 173 × 2,104 × 435 mm
Net weight	14.1 lbs / 6.4 kg
Measuring range	1 ft 12 in – 6 ft 7 in / 60 – 200 cm
Graduation	1/8 in / 1 mm

[Go to product website](#)



seca 453

Rechargeable battery pack for seca Scale-up Line*

- + Upgrade your Scale-up Line to the Mobile Version: compatible from serial number 10000001350445
- + For more than 8 hours of continuous weighing
- + Record weight and height in patient rooms, at the bedside, or wherever required
- + Fully charged within 3 hours

* Can only be retrofitted to seca Scale-up Line devices with serial number 10000001350445 or higher.

Dimensions (WxHxD)	3.7 × 13.8 × 2 in / 95 × 350 × 50 mm
Net weight	1.7 lbs / 0.8 kg
Cable length	29.5 in / 0.8 m

[Go to product website](#)

[Back to product overview](#)



seca 459

Panda bear Bao Bao for seca measuring stations and column scales

- + The cute little panda makes the measuring process more appealing and child friendly
- + Quick and easy installation
- + Easy to clean, with rounded edges
- + Long durability through high-quality DIBOND® material

Dimensions (WxHxD)	14.9 × 26.3 × 1.7 in / 380 × 670 × 44 mm
Net weight	1.1 lbs / 0.5 kg

[Go to product website](#)



seca 463

Barcode scanner mount

- + Ensures a barcode scanner is always ready at hand
- + Straightforward and quick mounting on the column
- + Suitable for all standard barcode scanners

* Barcode scanner from a third-party vendor

Dimensions (WxHxD)	3 × 3 × 6 in / 76 × 75 × 152 mm
Net weight	0.4 lbs / 0.2 kg

[Go to product website](#)



seca 487

Animal sticker for seca ultrasonic measuring stations

- + Colourful animals to make measuring stations more child-friendly
- + Animal figures at eye-level help obtain a correct head position for height measurement
- + High quality: sticker stays colorful for a long time
- + Approved medical device: Robust and resistant to disinfectants
- + Self-adhesive for easy application

Dimensions (WxHxD)	2.4 × 92.4 × 0 in / 60 × 2,346 × 0.5 mm
Net weight	0.4 lbs / 0.2 kg

[Go to product website](#)

[Back to product overview](#)



seca 223

Telescopic measuring rod for seca handrail scale

- + Extra-broad slide unit
- + Includes heel positioner for ideal three-point measurement
- + Developed for wall installation
- + Results are easy to read even for very tall patients
- + The rod is made of durable aluminum

Dimensions (WxHxD)	4.7 × 91.3 × 11.3 in / 120 × 2,320 × 286 mm
Net weight	2.1 lbs / 0.95 kg
Measuring range	2 – 90 in / 6 – 230 cm
Graduation	1/8 in / 1 mm

[Go to product website](#)



seca 452

Interface module to transmit measurements to EMR systems

- + Data transmission via WLAN or LAN
- + Patient identification right on device
- + No transcription errors

Dimensions (WxHxD)	3.7 × 4.4 × 1.4 in / 95 × 113 × 35 mm
Net weight	0.4 lbs / 0.18 kg

[Go to product website](#)

[Back to product overview](#)



seca 452

Interface module to transmit measurements to EMR systems

- + Data transmission via WLAN or LAN
- + Patient identification right on device
- + No transcription errors

Dimensions (WxHxD)	3.7 × 4.4 × 1.4 in / 95 × 113 × 35 mm
Net weight	0.4 lbs / 0.18 kg

[Go to product website](#)



seca 472

Stand for display elements of seca scales

- + Stable platform with non-slip rubber feet
- + Sturdy, durable construction
- + Can also be used as a stand for external printer

Dimensions (WxHxD)	18.1 × 40.7 × 18.1 in / 460 × 1,035 × 460 mm
Net weight	17.2 lbs / 7.8 kg

[Go to product website](#)



seca 405

Third ramp for seca XLine scales

- + Third ramp allows for easy entry from multiple sides
- + Sturdy build ensures long-lasting, accurate weigh-ins
- + Designed for hassle-free setup, ensuring quick usability
- + Features a yellow stripe marking the ramp's end for safety

Dimensions (WxHxD)	60 × 2 × 9.2 in / 1,524 × 52 × 234 mm
Net weight	17.4 lbs / 7.9 kg

[Back to product overview](#)



seca 452

Interface module to transmit measurements to EMR systems

- + Data transmission via WLAN or LAN
- + Patient identification right on device
- + No transcription errors

Dimensions (WxHxD)	3.7 × 4.4 × 1.4 in / 95 × 113 × 35 mm
Net weight	0.4 lbs / 0.18 kg



seca 454

Rechargeable battery for mobile use of EMR-integrated seca scales

Net weight	0.8 lbs / 0.35 kg
Cable length	45.2 in / 1.2 m

[Go to product website](#)

[Back to product overview](#)

Medical Measuring Systems and Scales since 1840

seca operates worldwide with
headquarters in Germany and
branches in:

seca france

seca united kingdom

seca north america

seca schweiz

seca zhong guo

seca nihon

seca mexico

seca austria

seca polska

seca middle east

seca suomi

seca américa latina

seca asia pacific

seca danmark

seca benelux

seca lietuva

and with exclusive partners in
more than 110 countries.

seca north america | west
13601 Benson Avenue
Chino, CA 91710
USA

seca north america | east
1439 Stuart Engals Blvd, Suite 200
Mount Pleasant, SC 29464
USA

T +1 800 542 7322
F +1 888 705 7397
E info.us@seca.com

Precision or health
[seca.com](https://www.seca.com)